

Modified protocol for Isolating RNA from Grape Stems

1. Ground sample with liquid nitrogen.
2. Weight 10-50 mg ground sample into a tube contains 400 μ l CPL Buffer. Mix the sample thoroughly by vortexing the tube 20 seconds.
3. (Optional) Incubate at 55-60C for 5 minutes.
4. Add 400 μ l chloroform and mix thoroughly by vortexing for 20 seconds.
5. Centrifuge at $\geq 13,000 \times g$ for 5 minutes.
6. Transfer 300 μ l of upper aqueous phase into a new 1.5 ml microtube.
7. Add equal volume of the RB Buffer followed by 325 μ l absolute ethanol (96-100%). Mix the sample thoroughly by vortexing.
8. Load 600 μ l of the sample into the HiBind RNA column. Centrifuge at 10,000 $\times g$ for 1 minutes. Discard the flow-through and re-use the collection tube.
9. Apply remaining sample to the HiBind RNA column by repeating step 8. Discard the flow-through and collection tube.
10. Place the RNA column into a new collection tube. Add 500 μ l RNA Wash Buffer I into the column . Centrifuge at 10,000 $\times g$ for 1 minutes. Discard the flow-through and re-use the collection tube.
11. Add 500 μ l RNA Wash Buffer II into the column . Centrifuge at 10,000 $\times g$ for 1 minutes. Discard the flow-through and re-use the collection tube.
12. Add 500 μ l RNA Wash Buffer II into the column . Centrifuge at 10,000 $\times g$ for 1 minutes. Discard the flow-through and re-use the collection tube.
13. Centrifuge the empty column for 2 minutes to dry the column.
14. Place the column into a new 1.5 ml centrifuge tube. Add 100 μ l DEPC Water into the center of the membrane in the column. Incubate at RT for 2 minutes.
15. Centrifuge at maximum speed ($>13,000 \times g$) for 1 minutes to elute RNA.

Modified protocol for Isolating RNA from Grape Leaves

1. Ground sample with liquid nitrogen.
2. Weight 10-50 mg ground sample into a tube contains 400 μ l CPL Buffer. Mix the sample thoroughly by vortexing the tube 20 seconds.
3. Add 400 μ l chloroform and mix thoroughly by vortexing for 20 seconds.
4. Centrifuge at $\geq 13,000 \times g$ for 5 minutes.
5. Transfer 300 μ l of upper aqueous phase into a new 1.5 ml microtube. Add 300 μ l of isoproanol to precipitate nucleic acid.
6. Centrifuge at 13,000 $\times g$ for 5 minutes. Remove the supernatant and dissolve RNA by 300 μ l DEPC Water.
7. Add equal volume of the RB Buffer followed by 325 μ l absolute ethanol (96-100%). Mix the sample thoroughly by vortexing.
8. Load 600 μ l of the sample into the HiBind RNA column. Centrifuge at 10,000 $\times g$ for 1 minutes. Discard the flow-through and re-use the collection tube.
9. Apply remaining sample to the HiBind RNA column by repeating step 8. Discard the flow-through and collection tube.
10. Place the RNA column into a new collection tube. Add 500 μ l RNA Wash Buffer I into the column . Centrifuge at 10,000 $\times g$ for 1 minutes. Discard the flow-through and re-use the collection tube.
11. Add 500 μ l RNA Wash Buffer II into the column . Centrifuge at 10,000 $\times g$ for 1 minutes. Discard the flow-through and re-use the collection tube.
12. Add 500 μ l RNA Wash Buffer II into the column . Centrifuge at 10,000 $\times g$ for 1 minutes. Discard the flow-through and re-use the collection tube.
13. Centrifuge the empty column for 2 minutes to dry the column.
14. Place the column into a new 1.5 ml centrifuge tube. Add 100 μ l DEPC Water into the center of the membrane in the column. Incubate at RT for 2 minutes.
15. Centrifuge at maximum speed ($>13,000 \times g$) for 1 minutes to elute RNA.