

## Contents

Introduction. . . . .	2
Benefits . . . . .	2
Binding Capacity. . . . .	2
Kit Contents. . . . .	3
Storage and Stability. . . . .	3
Before Starting. . . . .	3
E.Z.N.A.™ Cycle-Pure Kit Spin Protocol. . . . .	4
E.Z.N.A.™ Cycle-Pure Kit Vacuum/Spin Protocol. . . . .	5
Troubleshooting Guide. . . . .	6
Ordering Information. . . . .	7

## Introduction

The E.Z.N.A.™ family of products is an innovative system that radically simplifies the extraction and purification of nucleic acids from a variety of sources. The key to this system is the new HiBind® matrix that specifically, but reversibly, binds DNA or RNA under certain optimal conditions allowing proteins and other contaminants to be removed. Nucleic acids are easily eluted with deionized water or a low salt buffer.

The E.Z.N.A.™ Cycle-Pure Kit is a convenient system for the fast and reliable purification of PCR products. The E.Z.N.A.™ Cycle-Pure Kit uses HiBind® DNA technology to recover DNA bands from 100bp to 10kb free of oligonucleotides, nucleotides, and polymerase in yields exceeding 80%. Binding conditions are adjusted by the addition of a specially formulated buffer, and the sample is applied to a HiBind® DNA column. Following a rapid wash step, DNA is eluted with deionized water or a low salt buffer. Purified DNA is suitable for any downstream applications. No organic extractions or alcohol precipitations signifies a safe and rapid processing of multiple samples in parallel. Purified DNA can be directly used for most downstream applications include T-A ligations, PCR sequencing, restriction enzyme digestion, or various labeling reactions.

## Benefits

- **Fast**-DNA recovery from enzymatic reactions in less than 10 min
- **Reliability**-Optimized buffers that guarantee pure DNA every time
- **Safety**-No organic extractions
- **Quality**-Purified DNA is suitable for most applications

## Binding Capacity

Each HiBind® DNA column can bind ~30µg DNA.

## Kit Contents

Product Number	D6492-00	D6492-01	D649 2-02
Purification Times	5 preps	5 0 preps	200 preps
HiBind® DNA Columns	5	50	200
2 ml Collection Tubes	5	50	200
Buffer CP*	5 ml	40 ml	120 ml
Elution Buffer *	5 ml	10 ml	20 ml
DNA Wash Buffer*	1.5 ml	15 ml	3 x 25 ml
Instruction Booklet	1	1	1

\* Elution Buffer = 10mM Tris-HCl, pH 8.5

\* Buffer CP contains chaotropic salts which are irritants. The Equilibration Buffer contains Sodium Hydroxide. Take appropriate laboratory safety measures and wear gloves when handling.

\*The volume of the DNA Wash Buffer has been changed, see bottle label for dilution instruction.

## Storage and Stability

All E.Z.N.A.™ Cycle-Pure Kit components are guaranteed for at least 24 months from the date of purchase when stored at 22-25°C. Under cool ambient conditions crystals may form in Buffer CP. Simply warm to 37°C to dissolve.

## Before Starting

Please read the entire booklet to become familiar with the E.Z.N.A.™ Cycle Pure DNA Kit procedure

- DNA Wash Buffer must be diluted with **absolute ethanol** (96-100%) as follows and stored at room temperature.

D6492/D6493-00	Add 6ml of absolute ethanol to bottle
D6492/D6493-01	Add 60ml of absolute ethanol to bottle
D6492/D6493-02	Add 100ml of absolute ethanol to each bottle

## E.Z.N.A.™ Cycle-Pure Kit Spin Protocol

It is strongly advised that you familiarize yourself with the entire procedure before beginning this protocol. Omega Bio-Tek, Inc.'s E.Z.N.A.™ Cycle-Pure Kit is designed to be simple, fast, and reliable provided that all steps are followed diligently.

### Materials Supplied by User

- Microcentrifuge capable of at least 13,000x g
- Nuclease-free 1.5ml centrifuge tubes
- Optional: Sterile deionized water
- Absolute ethanol (~ 96-100%)

\*All centrifugation steps must be performed at room temperature

\*DNA Wash Buffer must be diluted with absolute ethanol prior to use (see page 3, or label instructions)

1. **Perform agarose gel/ethidium bromide electrophoresis to analyze PCR product.**
2. **Determine the volume of the PCR reaction, transfer the sample into a clean 1.5ml microcentrifuge tube, and add 4-5 volumes of Buffer CP.** For PCR products < 200bp add 6 volumes of Buffer CP.
3. **Vortex thoroughly to mix. Briefly spin the tube to collect any drops from the inside of the lid.**
4. **Place a HiBind® DNA column in a provided 2ml collection tube. Add 100µl Equilibration Buffer to the column. Incubate at room temperature for 4 minutes. Spin at maximum speed for 20 seconds.**
5. **Apply the sample to the HiBind® DNA column and centrifuge at 10,000x g for 1 min at room temperature. Discard the flow-through.**
6. **Wash the HiBind® DNA column by adding 700µl of DNA Wash Buffer diluted with absolute ethanol and centrifuge as above.**  
**IMPORTANT:** DNA Wash Buffer must be diluted with absolute ethanol before use. Refer to label for instructions. If refrigerated, DNA Wash Buffer must be brought to room temperature before use.
7. **Discard liquid and repeat Step 7 using 500µl of DNA Wash Buffer.**
8. **Discard liquid and centrifuge the empty HiBind® DNA column for 2 min at maximum speed (≥13,000x g) to dry the column matrix.** This is critical for good yields.
9. **Place HiBind® DNA column into a clean 1.5ml microcentrifuge tube. Add 30-50µl (depending on desired concentration of final product) of Elution Buffer (10mM Tris,**

**pH8.5) or water directly onto the column matrix and centrifuge for 1 min at  $\geq 13,000\times g$  to elute DNA.** This represents approximately 80-90% of bound DNA. An optional second elution will yield any residual DNA, though at a lower concentration.

- Yield and quality of DNA: Determine the absorbance of an appropriate dilution of the sample at 260nm and then at 280nm.** The DNA concentration is calculated as follows:

$$\text{DNA concentration} = A_{260} \times 50 \times (\text{Dilution Factor}) \mu\text{g/ml}$$

Fragments greater than 500bp in length can routinely be purified at  $> 80\%$  yield. Bands ranging from 100bp to 500bp gives yields of 60%- 90%. The ratio of ( $A_{260}$ ) / ( $A_{280}$ ) is an indication of nucleic acid purity. Alternatively, yield (as well as quality) can sometimes be best determined by agarose gel/ethidium bromide electrophoresis.

### E.Z.N.A.™ Cycle-Pure Vacuum/Spin Protocol

- Prepare the sample by following steps 1-3 of the spin protocol on page 4.**
- Prepare the vacuum manifold according to manufacturer's instructions.**
- Load the PCR reaction/CP solution from step 3 (of spin protocol) to the HiBind® DNA Column by decanting or pipetting and apply vacuum. After the samples have passed through the column switch off the vacuum source.**
- Wash the HiBind® DNA Column by adding 700µl of DNA Wash Buffer and turning on the vacuum source.**  
**IMPORTANT:** DNA Wash Buffer must be diluted with absolute ethanol before use. Refer to label for instructions.
- Repeat step 4.**
- Assemble the column into a 2ml collection tube and spin for 2 min at maximum speed ( $\geq 13,000\times g$ ) to dry the HiBind® DNA Column.**
- Place the column in a clean 1.5ml tube and add 30-50 µl of Elution Buffer(10mM Tris, pH8.5) directly onto the column matrix. Let it sit at room temperature for 1-2 minutes. Centrifuge for 1 min at  $\geq 13,000\times g$  to elute DNA.**

### Troubleshooting Guide

Problem	Solution
<b>Low DNA yields</b>	
Not enough Buffer CP added to sample	Add more Buffer CP as indicated. For DNA fragments $< 200\text{bp}$ in size, add up to 6 x volumes of Buffer CP.
Water pH is too low ( $< 7.5$ )	Check the pH of the water, adjust the pH of the water to 8.0 using Tris-HCl (2M, pH 8.5)
<b>No DNA eluted</b>	
DNA Wash Buffer has not been diluted with absolute ethanol (96-100%)	Prepare DNA Wash Buffer as instructed on the bottle, or refer to page 3.
<b>Optical densities do not agree with DNA yield on agarose gel</b>	
Trace contaminants were eluted from the column, thereby increasing $A_{260}$ .	Make Sure to Wash Column as instructed in steps 6 and 7 of the Spin Protocol, and steps 4 and 5 of the Vacuum/Spin Protocol. Alternatively, rely on agarose gel/ethidium bromide electrophoresis for quantization
<b>DNA sample floats out of well while loading agarose gel</b>	
Ethanol was not completely removed from column following wash steps	Centrifuge column as instructed in step 8 of the spin protocol and step 6 of the Vacuum/Spin protocol.

### Ordering Information

Product	Applications	Cat. No.
Cycle-Pure Kit	PCR product purification	D6493-01/02 D6492-01/02
MicroElute™ Cycle-Pure Kit	PCR product purification - special column for lower elution volume	D6293-01/02
E-Z 96® Cycle-Pure Kit	PCR product purification in a 96-well format	D1043-01/02
Mag-Bind® Cycle-Pure Kit	PCR product purification with magnetic beads	M1322-01/02
MicroElute™ DNA Cleanup Kit	DNA recovery from enzymatic reactions- special column for lower elution volume	D6296-01/02
Mag-Bind® Oligonucleotide Purification Kit	DNA recovery from enzymatic reactions using magnetic beads	M2514-01/02
E-Z 96® Mag-Bind® Sequencing Dye Removal Kit	Sequencing dye terminator removal with magnetic beads in a 96-well format	M1320-01/02
DNA Probe Purification Kit	DNA Cleanup from enzymatic reactions	D6538-01/02
Gel Extraction Kit	DNA recovery from Agarose Gel	D2501-01/02 D2500-01/02
MicroElute™ Gel Extraction Kit	DNA recovery from agarose gel- special column for lower elution volume	D6294-01/02
Ultra-Sep® Gel Extraction Kit	Low cost DNA recovery from agarose gel	D2510-01/02
Poly-Gel DNA Extraction Kit	DNA recovery from poly-acrylamide gel	D2561-01/02
Poly-Gel RNA Extraction Kit	RNA recovery from poly-acrylamide gel	R6376-01/02
RNA Probe Purification Kit	RNA Cleanup from enzymatic reactions	R6249-01/02

MicroElute™ RNA Cleanup Kit	RNA recovery from enzymatic reactions- special column for lower elution volume	R6247-01/02
-----------------------------	--	-------------

Product	Size	Product No.
Buffer CP	200ml/500ml	PDR042/PDR043
Elution Buffer	100 ml	PDR048
DNA Wash Buffer	(40 ml; add 60ml ETOH)	PDR044
DNA Wash Buffer	(100 ml; add 150ml ETOH)	PS010
DNA Wash Buffer	(500 ml; add 750 ETOH)	PS011
2ml capless collection tubes	500/BAG	SS1-1370-00
1.5ml DNase/RNase Free Centrifuge Tubes	500/BAG	SS1-1210-00

Please Call, Fax , or e-mail us to place an order.

Tel: 770-931-8400 (US) Fax: 770-931-0230 (US) e-mail: [info@omegabiotek.com](mailto:info@omegabiotek.com)

Tel: 1-800-832-8896 (Toll free) Fax: 1-888-624-1688 (Toll free)

[www.omegabiotek.com](http://www.omegabiotek.com)